

SoloSpecTM Sigmoidoscope

Orion Concepts Inc's disposable Sigmoidoscopes are manufactured in a clear environmentally friendly plastic with excellent light transmission which gives an optimal, distal, annular light view for the practitioner. The Sigmoidoscope has smooth edges that contribute to a gentle examination for the patient.

Quality: ISO 13485


Standard size*

Small size*



Benefits include

- Contoured finger grip for ease of use
- Clearly marked graduated scale is easily read and provides a smooth surface
- Fits most popular fiber-optic light-heads/illuminators
- A smooth transition between the tip of the obturator and the end of the Sigmoidoscope ensures maximum patient comfort
- Exceptional distal tip illumination gives a bright and clear view



Additional information

- Single-use only
- Not autoclaveable
- Environmentally friendly

Order Code:

Small size: 42319-1C* Standard size: 42319-2C*

*Sold in case of 25

Distributed by:



orionconceptsinc.com
info@orionconceptsinc.com