

Sureband[™]

Introducing the multi-band ligator


that's latex-free, simple, fast, reliable
and ready to use.



More clinically efficient and easier to use, Sureband™ is better for both surgeons and patients.

Having been designed by surgeons for surgeons, the Sureband™ makes ligation procedures more comfortable for patients, more clinically efficient for surgeons and easier for nurses.

Designed for use with an
ordinary anoscope



Innovative multi-action
handle controls suction,
fires and reloads bands

Clinically more efficient

Sureband™ is very clinically efficient. Its multi-action handle means you only need one hand to control suction and release and reload bands. This leaves your other hand free to hold a proctoscope and means you can work alone. The design of Sureband™ also allows for much more accurate applications, reducing the need for additional procedures.

Better patient outcomes

With the shortened procedure times achieved and latex-free bands, the Sureband™ improves patient comfort considerably and its disposability eliminates the risk of cross-contamination.

Ready to use

The Sureband™ gun is ready to use straight from the airtight, clean tray. Kitted with an anoscope complete with light guide, there are no spare parts to lose, find or fit and the only extra equipment you need is standard hospital suction and an illuminator.

Easier and faster

Because the gun is pre-loaded, already clean and can be used single-handed, the time taken for a typical ligation procedure can be reduced from 20 minutes to between 2 and 4 minutes for the banding of up to 3 haemorrhoids. In addition, if purchased in kit form, the illuminated anoscope delivers light to exactly where it's needed.



Pre-loaded with 4 latex-free bands. The deep nozzle offers optimum suction pressure for more accurate banding

Slim 120mm insertion probe

Sureband
ORION

The Revolutionary Multi-Band Ligator

available with illuminated anoscope

A Transparent Nozzle

The clarity of the tip maximizes visibility throughout the entire procedure. Preloaded with 4 green, latex-free, bands. The deep nozzle permits optimum suction pressure for more accurate banding.

Latex-Free Green Molded Bands

Each band is individually molded to create round edges meaning that the band rolls off easily without twisting and ensures a more reliable release & secure fit onto the tissue.

Innovative multi-action handle controls suction, fires & reloads bands. The stepped trigger line guides the user through each step of the application.

Sureband™ (Ligator only)

Part No : 40810 - Case Qty/10

Sureband™ Kit

Includes Sureband™ and disposable medium anoscope with light guide*

Part No : 40811C - Case Qty/5


Sureband™ Kit

Includes Sureband™ and disposable large anoscope with light guide*

Part No : 40811-3C - Case Qty/5

* LED illumination handle (GS-LEI-A-HB) sold separately

Distributed by:



Slim 120mm barrel



orionconceptsinc.com

Manufacturers of GreenStar™ Endoscopic Illuminator