


## SoloSpec<sup>™</sup> Anoscope with Bevel tip

Introducing the new SoloSpec<sup>™</sup>  
disposable anoscope with bevel tip.


Made of environmentally friendly plastic,  
it conducts light very well and provides  
an optimal view for the examiner and  
a gentle examination for the patient.





## Benefits

- The beveled tip provides better access for treatment of hemorrhoids.
- Smooth transition between the obturator and anoscope body
- Compatible with Welch Allyn & Heine Illumination systems
- Improves infection control by reducing potential patient to patient cross-contamination



## Material

- Polystyrene crystal, Polyethylene HD
- Environmentally friendly

### Order Code:

42310-2-SC – SoloSpec™ Disposable Anoscope Bevel Tip. Case /25

Distributed by:

**ORION**<sup>™</sup>  
CONCEPTS

orionconceptsinc.com  
info@orionconceptsinc.com